

**Doug Juers**, Sean Stothers,  
Jack Quire, Charlie Dibb  
Department of Physics & Program in BBMB  
Whitman College, Walla Walla WA


## Background

**Macromolecular crystals** are easily damaged during transfer to new solutions due to handling and/or osmotic effects. Here we explore two methods to gently equilibrate crystals to new solutions with minimal handling.


## Approach

old solution  
new solution

## Liquid Based



## Aerosol Based





## References/Acknowledgments

1. Sina Boeshaghi, Eduardo da Veiga Beltrame, Dylan Bannon, Jase Gehring and Lior Pachter, [Principles of open source bioinstrumentation applied to the poseidon syringe pump system](#), Scientific Reports 9, Article number: 12385 (2019)
2. This work was supported in part by a grant from the National Institutes of Health (GM0900248).



# Gentle equilibration of macromolecular crystals to new solutions with minimal handling.

## Results


### Liquid Based




Top: 40 µL pot (pcr tube cap) with inflow and outflow tubes. Bottom: view from below showing crystal on right.




Top: system showing open-source syringe pumping system<sup>1</sup> with enhancements for crystal equilibration. Bottom: auto gradient view of control software



Alpha lactalbumin  
0 -> 25 % glycerol  
Direct (left) vs 50 min linear gradient




Thermolysin  
~2 M AmSO<sub>4</sub> -> 0 M Am SO<sub>4</sub>  
Direct (left) vs 15 min linear gradient



Lysozyme  
8 % NaCl -> 3% NaCl  
Direct (left 2 xtals) vs 40 min linear gradient

## Aerosol Based

### Batch/Closed System Method




Draw vape-pen generated aerosol from syringe





Eject aerosol into cryovial with loaded sample and incubate.


### Process/Open System Method



Cryoprotectant is directed towards a loop mounted crystals sitting in a cryovial with a nebulizer plate.



Thermolysin  
~2 M AmSO<sub>4</sub> -> 0 M Am SO<sub>4</sub>  
Direct (left) vs ~10 min incubation in vial



Alpha lactalbumin  
0 -> 25 % glycerol  
Direct (left) vs ~5 min incubation in vial

## Summary

**The liquid based approach** has proven to be robust and reliable.

**The aerosol based approach** has been shown to work in principle. The many process variables make it less reliable. However it has potential to be a more rapid method, and further development is underway.